
# Liquid Nitrogen(LN2) Automatic Supply Device (Demas)

Manual ver.1

**Safety first !**

①

Put your container  
on the weight scale.

You can use this device only if  
ordered on Cryogen Ordering  
System.

②

Check if the screen is  
the same as below.

If the screen saver is  
displayed, touch it once.

③

Scan the barcode on  
your container with a  
barcode reader.

!

If you place multiple orders, please choose the appropriate one.

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④

Press "Weighing Button ●" and measure your container weight.

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AIR LIQUIDE JAPAN

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⑤

Insert the filling nozzle into your container and press the "Confirm Button".

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AIR LIQUIDE JAPAN

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⑥

Press "Full Filling Button(1)" or "Volume Setting Button(2)" to set the required amount.

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AIR LIQUIDE JAPAN

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\* Max filling capacity is 90% of container capacity.

⑦

Press "Filling Start Button ●" and start LN2 filling.

⑧


Until the piping is completely cooled, release nitrogen gas from the discharge valve.

Gas is released from the piping behind. Please wait for the pipes to cool down.

⑨

When the piping has cooled down, it switches to automatic filling and starts filling the container with LN2.

Wait for LN2 filling to complete.

!

LN2 filling amount is measured with a weight scale.

Do not step on the scale or place objects during filling.

When it is detected that the weight is larger than the set value, LN2 filling is completed. However, since LN2 filling amount is determined on step 11, LN2 filling amount does not increase.

⑩


When LN2 filling is completed, an alarm sounds and press "Alarm Stop Button".

Even if LN2 filling is completed, do not move the container yet.

11

Remove the filling nozzle from the container and press "Weighing Button ●".

The screenshot shows the DeMaS (液化窒素自動充填システム) interface. At the top, it displays the system name and the time 13:45:38. Below this, it shows the tare weight (ロードセル重量) as 7.0 kg. The user information (利用者) is listed as test 上田 (氏名) and 低温センター (所属). The elapsed time (経過時間) is 3分21秒. The main message states "充填終了" (Filling completed) and "充填作業が終了しました。充填量は 4.1Lです。" (Filling operation has ended. The filling amount is 4.1L). At the bottom, there is an AIR LIQUIDE JAPAN logo.

The final filling weight is determined by the weight at the time of pressing "Weighing Button". This amount is automatically entered into the Cryogen Ordering System.

12

LN2 filling is completed.

Within 30 minutes after filling, it can be refilled. If refilling, start from step 1.

If you have any questions, please inform the staff.